


The background features a light blue-grey grid pattern. On the left side, there are several overlapping geometric shapes, including triangles and circles, in shades of yellow, orange, and grey. On the right side, there are concentric circles in shades of orange and grey. The overall aesthetic is modern and technical.

Joint Angles Visualization in Running Exercise

1.1. Phases of Sprinting



Which Muscle Groups Are Used in Running?

How the Hip Moves while Running

1. Introduction

1.2. Recording System

1.2.1. Pressure-Instrumented Treadmills - myoPRESSURE™


Individually-calibrated capacitive sensors measure data to conduct static and dynamic plantar pressure mapping.



1. Introduction

1.2. Recording System

1.2.1. Pressure-Instrumented Treadmills - myoPRESSURE™



AutoSave Off

File Home Insert Page Layout Formulas Data

Paste

Clipboard Font

N16

MR32023-07-21-16-30

Ready Accessibility Unavailable

1	A	B	C	D	E	F	G	H	I	J	K	L
2	type	begin_time	frequency	count	created wi	exported w	project	last_name	first_name	born	sex	measurement_date
3	record	0	100	8560	Noraxon M3.20.18	operator te	Hoang Trung Hieu				male	2023-07-21T09:30:08.334Z
4	time	Activity	Marker	LT Max Pre	RT Max Pressure (N/cm?)					Trials in Ne		
5	0	Activity		0	0							
6	0.01			0	0							
7	0.02			0	0							
8	0.03			0	0							
9	0.04			0	0							
10	0.05			0	0							
11	0.06			0	0							
12	0.07			0	0							
13	0.08			0	0							
14	0.09			0	0							

Frame

Left and Right Max
pressure on each
foot

Start time of each
cutted file



Start time of this
file

The two legs should alternate with each other like this

1. Introduction

1.2. Recording System

1.2.2. Motion Capture for Biomechanics - Nexus



3D coordinates of
markers

1.2.2. Motion Capture for Biomechanics - Nexus

NexusHieu running04

NexusHieu running05

NexusHieu running06

NexusHieu running07

NexusHieu running08

	A	B	C	D	E	F	G	H
1								
2	Model Outputs							
3	100							
4	Frame	Sub Frame X	Hoang Trung Hieu lower pig:LAbsAnkleAngle	Y	Z	Hoang Trung Hieu lower pig:LAnkleAngles		
5		deg		deg	deg	X	Y	Z
6						deg	deg	deg
7	16	0	-3.99616	0	0	-3.99639	1.37895	-0.603
8	17	0	-5.37514	0	0	-5.38504	1.78641	-3.469
9	18	0	-6.55811	0	0	-6.57601	1.89366	-4.221
10	19	0	-8.65248	0	0	-8.66401	1.71238	-2.949
11	20	0	-11.0276	0	0	-11.0276	1.3057	-0.086
12	21	0	-13.1916	0	0	-13.2175	0.788645	3.551
13	22	0	-15.4283	0	0	-15.5502	0.281637	7.091
14	23	0	-17.7914	0	0	-18.0816	-0.15692	10.11



Frame

Angle name
(L: Left + Angle angle)

Angle in


1. Introduction

1.2. Recording System




Simulation of file segments of Nexus corresponding to MR file in Frame field t_1, t_2, \dots : start frame of n -th segment of Nexus file

2. Workflow




2.1. Process (1): Processing the MR file and receive the information of each phase in cycles

- First, we need to assign the process to the stick bar for easier visualization.
- The file is distributed longitudinally, so the longitudinal stick bar represents the pressure values that differ from 0 in two columns: Left foot pressure and Right foot pressure.
- The gap between each bar in the longitudinal direction represents values of 0.
- The frame increases from the top to the bottom of the chart.
- The main task of these bars is to identify the time when the pressure changes from 0 to a non-zero value. This is the moment when the phase changes.



Function 1: Find value come back from pressure to zero : find_zero_seg



Input :



- pressure : array of column of pressure
- frame_start: first frame of pressure sequence

Output:

- The last frame of pressure sequence

This function aims to find what frame that the pressure sequence come back to zero sequence. This point can be end of stance phase or end of mid swing phase

frame_end_pressure



Function 2: Find value come back from zero to pressure: find_pres_seg

Input :

- pressure : array of column of pressure
- frame_end_pressure: last frame of pressure sequence


Output:

- The last frame of zero sequence


*The structure is inverted
with Function 1*

This function aims to find what frame that the zero sequence come back to pressure sequence. This point can be end of stance phase or end of mid swing phase

frame_end_pressure





Function 3: Receive the phase information from MR file



- This function aims to determine which foot starts first.
- If the Left foot starts first, it means the value in the left pressure column differs from 0 first. The flag `check_start_left_or_right` is set to 0, and then Block A is executed.
- Block A loops through each row in the Excel file, returning a list containing the stance, initial swing, mid swing, and terminal swing phases of each foot, and saves the results to another Excel file.
- Otherwise, if the Right foot starts first, `check_start_left_or_right` is set to 1, and Block B is executed.
- The difference between Block A and Block B is that Block A is designed for cases where the Left foot starts first, while Block B is designed for cases where the Right foot starts first.

Block A




2.2. Process (2) and (3): Find index of angle on Nexus file with corresponding frame on MR file

Define the initial construct

First, we need define the corressponding angle on each frame on Nexus file with the corresponding pressure on each frame on MR file. Because MR file contrain information about the cycles

Let imagine that initial we have the whole MR file, so we represent the frame column of MR file by a long green stick

Nexus produce several files, we represent it by several shot stick, with blue stick is frame column of Nexus file, yellow stick is corressponding angle at that frame



At the first time , there are a gap on 2 kind of file because we can not locate what frame of Nexus corresponding with what frame of MR

2.2. Process (2) and (3): Find index of angle on Nexus file with corresponding frame on MR file

Process (2):

Let read the start time of the each recording.

Assume


- t_0 is start time of MR file on second
- t_n is start time of Nexus segment file on second
- 1 frame = 1/100 second

Then the first frame of Nexus file will be attached with the frame $(t_n - t_0) * 100$ of MR file

Process (3):


Only keep the corressponding frame on MR with angle on Nexus

Can be considered as a DataFrame contains 3 column with NaN at the gap of 2nd and 3rd column



Can be considered as a DataFrame contains 2 column without NaN, 1st column is frame, 2nd column is angle

2.3. Process (4): Creating DataFrame contain angle on cycles





Start frame
of cycle i-th

End frame
of cycle i-th

Frame	Angle
...	...
...	...
...	...

DataFrame A





angles_in_cycles then have the form as the dataframe above, with each column represent each cycle and each row represent each angle at each frame on that cycle. So It often contain the NaN value at several bottom row of some column because the length of each cycle are usually difference

2.4. Process (5) and (6): Plotting chart from the calculated information

Process (5):

Plot the plot with x-axis represent the frame in cycle, y-axis represent the angle value at that frame. Each line is drawn from each column in DataFrame B this plot help to see the fluctuation of angle of person on running.


However, DataFrame B originally have the cycles with the difference length, we need to normalize the length to one fixed value to see the all the average representations of angles in a cycle. Then the length will be normalized from 0 to 100% and the fluctuation will be visualized by standard deviation of these point

Process (6): with each x from 1 to 100 (integer), we will take the angle at $x\%$ of each cycle, to calculate mean and std at point $x\%$

Assume:

- T_i is the i -th cycle are considering
- Iterate through each value x in $\text{range}(1,100)$
- Iterate through each column in DataFrame B, determine the frame in cycle which is at $x\%$ of cycle length T_i , extract the angle at this frame, append to list containing the set of angle at $x\%$
- The result is expected that at each $x\%$ from 1 to 100 will be contain a set of angle value

Here is an example of what happened in the process (6), with T_i (have length a), T_j (have length b) is two random cycle, they have difference value of cycle length although both lies on an DataFrame B, so the shorter cycle (T_j) will contain some NaN values. So $x\%$ of T_j and $x\%$ of T_i will have difference index on the frame column. This is why we need to re-calculate the frame index in the term of percent



3. Results:

